


Alkaline Phosphatase

End Point Determination

Store at 2- 8°C

PRINCIPLE

Colorimetric determination of the ALP activity according to the following reaction:



The phenol liberated is measured in the presence of amin 4 antipyrine and potassium ferricyanide.

The presence of sodium arsenate in the reagent stops the enzymatic reaction

REFERENCE VALUES

Adults	21 – 92 U/l
Children	71 – 142 U/l

Children have a higher normal value. It is strongly suggested that each laboratory establish its own normal range.

SAMPLES

Serum or heparinized plasma can be used. Stable for 7 days at 2-8°C. The activity increases if samples are left at room temperature for several hours. Hemolysis will interfere.

REAGENTS

R1:

Disodium phenylphosphate	5.0 mmol/l
Carbonate-bicarbonate buffer	50 mmol/l

R2:

Phenol	142 U/l
--------	---------

R3:

Amino 4 antipyrine	60 mmol/l
Sodium Arsenate	240 mmol/l
TOXIC. CORROSIVE REAGENT	
Use an automatic pipette.	
Avoid contact with the skin	

R4:

Potassium Ferricyanide	150 mmol/l
------------------------	------------

PREPARATION OF WORKING REAGENT

All reagents are ready to use.

PROCEDURE

Wavelength	510 nm (Hg 492)
Zero adjustment	Reagent Blank

Set up the following tubes:

	Serum sample	Serum Blank	Standard	Reagent Blank
Reagent 1	2 ml	2 ml	2 ml	2 ml

Incubation for 5 minutes at 37 C

Serum	50 µl	-	-	-
Reagent 2	-	-	50 µl	-

Incubation for exactly 15 minutes at 37 C

Reagent 3	0.5 ml	0.5 ml	0.5 ml	0.5 ml
-----------	--------	--------	--------	--------

Mix well

Reagent 4	0.5 ml	0.5 ml	0.5 ml	0.5 ml
-----------	--------	--------	--------	--------

Serum	-	50 µl	-	-
-------	---	-------	---	---

D-water	-	-	-	50 µl
---------	---	---	---	-------

Mix. Let stand for 10 minutes in the dark

Measure.

The color intensity is stable 45 minutes

CALCULATION

$$\text{ALP U/l} = \frac{(\text{A}) \text{ serum sample} - (\text{A}) \text{ serum Blank}}{(\text{A}) \text{ Standard}} \times 142$$

LINEARITY


Any sample that gives a value >285 U/L should be diluted.

BIBLIOGRAPHY

Kind P.R.N, King E.G J.Clin. Path 1954-7-322
BELFIELD A, GOLDBERG D.M Enzyme 1971-12-561

The following symbols are used on labels

 For in vitro diagnostic use

 Use day (last day of the month)

 Temperature limitation

 Batch code

 Code